

PHANTOM SERIES


BB STYLE CERVICAL RETRACTOR BLADES

Features

- Non-glare matte finish
- Wide range of blade lengths and widths
- Unique blade shape prevents climbing out of surgical site
- Multi-level retraction capability
- Titanium radiolucent



**TeDan Surgical
Innovations**



TSI BB Style Retractor Blades are compatible with all TeDan and competitor retractors that accept a ball snap blade.



CS-1320 through CS-1327 Cervical BB Style Blades with Long Teeth

- 22mm wide
- Available lengths: 34mm, 40mm, 45mm, 50mm, 55mm, 60mm, 70mm



CS-1330 through CS-1337 Cervical BB Style Blunt Blades

- 22mm wide
- Available lengths: 34mm, 40mm, 45mm, 50mm, 55mm, 60mm, 70mm



CS-1341 through CS-1348 Cervical BB Style Blades for Multi Level with Long Teeth

- 27mm wide
- Available lengths: 40mm, 47mm, 55mm, 65mm, 70mm, 75mm




CS-1351 through CS-1358 Cervical BB Style Blunt Blades for Multi Level

- 27mm wide
- Available lengths: 40mm, 47mm, 55mm, 65mm, 70mm, 75mm

Phantom Series

BB Style Cervical Retractor Blades



CS-1900

CS-1900 PHANTOM CS BB STYLE CERVICAL RETRACTOR SET

	ea.		ea.
CS-0100 Phantom CS Transverse Retractor with Double Hinge, 120mm rack	1	CS-1344 Phantom CS BB Style Cervical Blade Multi Level Long Teeth 27 x 55mm, Titanium	2
CS-0102 Phantom CS Transverse Step Over Retractor with Double Hinge, 120mm rack	1	CS-1346 Phantom CS BB Style Cervical Blade Multi Level Long Teeth 27 x 65mm, Titanium	2
CS-0105 Phantom CS Side-Loading Retractor Blade Handle	2	CS-1347 Phantom CS BB Style Cervical Blade Multi Level Long Teeth 27 x 70mm, Titanium	2
CS-1320 Phantom CS BB Style Cervical Blade Long Teeth 22 x 34mm, Titanium	2	CS-1348 Phantom CS BB Style Cervical Blade Multi Level Long Teeth 27 x 75mm, Titanium	2
CS-1321 Phantom CS BB Style Cervical Blade Long Teeth 22 x 40mm, Titanium	2	CS-1351 Phantom CS BB Style Cervical Blade Multi Level Blunt 27 x 40mm, Titanium	2
CS-1322 Phantom CS BB Style Cervical Blade Long Teeth 22 x 45mm, Titanium	2	CS-1352 Phantom CS BB Style Cervical Blade Multi Level Blunt 27 x 47mm, Titanium	2
CS-1323 Phantom CS BB Style Cervical Blade Long Teeth 22 x 50mm, Titanium	2	CS-1354 Phantom CS BB Style Cervical Blade Multi Level Blunt 27 x 55mm, Titanium	2
CS-1324 Phantom CS BB Style Cervical Blade Long Teeth 22 x 55mm, Titanium	2	CS-1356 Phantom CS BB Style Cervical Blade Multi Level Blunt 27 x 65mm, Titanium	2
CS-1325 Phantom CS BB Style Cervical Blade Long Teeth 22 x 60mm, Titanium	2	CS-1357 Phantom CS BB Style Cervical Blade Multi Level Blunt 27 x 70mm, Titanium	2
CS-1327 Phantom CS BB Style Cervical Blade Long Teeth 22 x 70mm, Titanium	2	CS-1358 Phantom CS BB Style Cervical Blade Multi Level Blunt 27 x 75mm, Titanium	2
CS-1330 Phantom CS BB Style Cervical Blade Blunt 22 x 34mm, Titanium	2	CS-0900 Phantom CS Sterilization Tray (for CS-1900 Set)	1
CS-1331 Phantom CS BB Style Cervical Blade Blunt 22 x 40mm, Titanium	2		
CS-1332 Phantom CS BB Style Cervical Blade Blunt 22 x 45mm, Titanium	2		
CS-1333 Phantom CS BB Style Cervical Blade Blunt 22 x 50mm, Titanium	2		
CS-1334 Phantom CS BB Style Cervical Blade Blunt 22 x 55mm, Titanium	2		
CS-1335 Phantom CS BB Style Cervical Blade Blunt 22 x 60mm, Titanium	2		
CS-1337 Phantom CS BB Style Cervical Blade Blunt 22 x 70mm, Titanium	2		
CS-1341 Phantom CS BB Style Cervical Blade Multi Level Long Teeth 27 x 40mm, Titanium	2		
CS-1342 Phantom CS BB Style Cervical Blade Multi Level Long Teeth 27 x 47mm, Titanium	2		

CONTACT US

TeDan Surgical Innovations, LLC
Ph: 713-726-0886 • Fax: 713-726-0846 • Toll Free: 877-726-0886
© TSI 2012

12615 W. Airport Blvd, Suite 200 • Sugar Land, TX 77478
www.tedansurgical.com • sales@tedansurgical.com

