


# DISTRACTION SCREWS FOR ALL ANATOMICAL NEEDS



- **Bullet Shaped Head**  
Will not damage the Holding Spring
- **Long Self-Tapping Flute**
- **Fits all Distraction Screw Drivers**
- **Available in Sterile and Non-Sterile Packaging**
- **Sharp Tip**  
Provides easy bone penetration


The Bullet shaped head fits in all distraction screw drivers without damaging the holding wire. A sharp tip with Self-Tapping threads allows for easy, free-hand placement. Available in 12, 14 & 16mm lengths in either sterile or non-sterile packaging.




Phantom CS™ Dual Drill Guides offer unique parallel distraction screw placement every time!



CS-0110  
Dual Distraction Screw Drill Guide



Foot Plate on Dual & Single Drill Guides align the placement of the guide tubes.



Distraction Screws parallel alignment guarantees parallel distraction without distractor slippage.

Quick and easy to penetrate the cortical bone for distraction screw placement.

Phantom CS™ Twist Drill  
CS-0109

Phantom CS™ Twist Drill with handle for quick easy drilling.



Drill and place distraction screws directly through the Drill Guides.

## ORDER INFORMATION

| Description | Sterile Product * | Non-Sterile ** |
|-------------|-------------------|----------------|
| 12mm        | DS-0012           | DS-0212        |
| 14mm        | DS-0014           | DS-0214        |
| 16mm        | DS-0016           | DS-0216        |

\* Each box contains 5 individually packaged units.  
\*\* One unit per bag.

## MANUFACTURED FOR

TeDan Surgical Innovations  
11333 Chimney Rock Road  
Suite #180  
Houston, Texas 77035

Ph: 713-726-0886 • Fax: 713-726-0846  
Toll Free: 877-726-0886  
www.tedansurgical.com  
Email: sales@tedansurgical.com