


# PHANTOM LS™ MINI LUMBAR RETRACTOR SYSTEM

**TSI** TeDan Surgical  
Innovations



LS-1704

**Phantom LS Mini Lumbar Retractor Set** facilitates retraction to access the posterior lumbar spine. Retractor blades are semi-radiolucent and ideal for procedures using fluoroscopy.



LS-0100 retractor arms allow for the attachment of multiple retractor blades

- 1 Adjustable retractor blade holder creates an elliptical incision
- 2 Unique extended length blade holder improves access to the pedicles
- 3 Hinged arms conform to patient's anatomy

## Phantom LS™ Mini Lumbar Retractor Set

facilitates retraction to access the posterior lumbar spine. Retractor blades are semi-radiolucent and ideal for procedures using fluoroscopy.

CAT. NO.	DESCRIPTION	QTY.
CS-0105	PHANTOM CS SIDE-LOADING RETRACTOR BLADE HANDLE	2
LS-0030	PHANTOM LS LUMBAR BLADE 15° W/TEETH, 23X35MM, TITANIUM	2
LS-0031	PHANTOM LS LUMBAR BLADE 15° W/TEETH, 23X45MM, TITANIUM	2
LS-0032	PHANTOM LS LUMBAR BLADE 15° W/TEETH, 23X55MM, TITANIUM	2
LS-0033	PHANTOM LS LUMBAR BLADE 15° W/TEETH, 23X65MM, TITANIUM	2
LS-0034	PHANTOM LS LUMBAR BLADE 15° W/TEETH, 23X75MM, TITANIUM	2
LS-0100	PHANTOM LS LUMBAR RETRACTOR BODY	1
LS-0101	PHANTOM LS T-HANDLE CRANK	1
LS-0102	PHANTOM LS ADJUSTABLE BLADE HOLDER	4
LS-0201	PHANTOM LT TOP LOADING RETRACTOR BLADE HOLDER	6
LS-0704	PHANTOM LS STERILIZATION TRAY (FOR LS-1704 MINI SETS)	1



## CONTACT US

TeDan Surgical Innovations, LLC  
12615 W. Airport Blvd, Suite 200 • Sugar Land, TX 77478  
Ph: 713-726-0886 • Fax: 713-726-0846 • Toll Free: 877-726-0886  
[www.tedansurgical.com](http://www.tedansurgical.com) • [sales@tedansurgical.com](mailto:sales@tedansurgical.com)

